

Oscillometry, the missing piece in lung transplant rejection

Transplanted lungs are at risk of rejection after transplant. Frequent surveillance, and early diagnosis are essential for clinical outcomes. Rejection is diagnosed primarily using spirometry and bronchial biopsy, but there is a need for additional, more sensitive, and less-invasive tests.

Oscillometry (OS) is a technology that measures lung mechanics during tidal breathing and has been shown recently to be a useful adjunct method to detect and monitor acute and chronic lung transplant rejection as well as other airway diseases¹.

Key measurements using OS:

Key Measurement	Description
R5	Resistance to airflow of the whole lung.
R5-19	Resistance heterogeneity of the airways. <i>Increased R5-19 signals potential small airway dysfunction.</i>
X5	Elastance of the respiratory system. <i>Decreased X5 implies increased tissue elastance.</i>
AX	Increased AX correlates with airspace de-recruitment, ventilation defects, and increased lung stiffness.

Acute cellular rejection (ACR)

ACR begins in the small airways (i.e. the 'silent zone'²). Currently, bronchial biopsy is the gold standard to diagnose ACR in addition to monitoring with spirometry at a patient's home or at pulmonary function labs.

Challenges present with conventional diagnostic methods:

Spirometry	<p>Spirometry lacks sensitivity to detect biopsy-proven episodes of ACR³.</p> <ul style="list-style-type: none">• Up to 75% of all the small airways must change before being detected by forced expiratory volume (FEV₁)⁴• Spirometry has 60% sensitivity for detecting clinically significant grade 2 ACR (10% drop in FEV₁) or higher⁵• Up to 25% of patients are asymptomatic and would only be diagnosed if surveillance bronchoscopy is performed^{6,7}
Bronchial Biopsy	<ul style="list-style-type: none">• Transbronchial biopsy is limited by the sampling of a small region of the lung and usually contains only a few small airways• Bronchoscopy and transbronchial biopsies have associated risks, procedure-related morbidity, and costs• Routine surveillance bronchoscopy with transbronchial biopsies is not universally adopted by lung transplant centers

Find the missing piece for ACR

OS has been shown to assess lung function including the small airways. It has the potential to detect subtle physiological changes associated with rejection in transplant patients, prior to symptomatic presentation or spirometry findings.

In a recent publication by Cho et al.³, it was shown that OS results can assist:



Early diagnosis - Abnormal OS parameters are detected 2 weeks prior to scheduled transbronchial biopsy. This is as early as 4 weeks post-transplant for patients with grade 2 ACR.



Surveillance - OS detects significant airway obstruction in double lung transplant patients with grade 2 ACR, while FEV₁ measurements were stable or improved for patients post-transplantation.



Treatment evaluation - OS measurements return to normal as early as 1 week after treatment starts, suggesting treatment efficacy can be monitored using OS.

"The data of Cho and colleagues support... the utility of OS as an essential adjunct to spirometry in [ACR] patients, particularly in centers that do not perform routine surveillance bronchoscopy with transbronchial biopsy."²

-Usmani et al. 2020

Chronic lung allograft dysfunction (CLAD)

Chronic rejection, or CLAD, develops in 50–70% of patients within 5 years and is the primary cause of death in recipients⁸. CLAD is diagnosed by a 20% decline in FEV₁ when compared to the average of the two of the best post-transplant measurements taken at least 3 weeks apart⁷.

Find the missing piece for CLAD



Patient Surveillance - Intra-subject variance of OS parameters (R5-19, X5, AX), but not FEV₁, was strongly and independently associated with an increased risk of CLAD in double lung transplant recipients measured a year post transplant⁸.



Diagnosis - In patients with high and low bronchial biopsy A-scores with similar FEV₁, intra-subject variance showed significant difference between the two groups, suggesting OS may be helpful in detecting CLAD in patients prior to the drop of FEV₁.

Additional research is ongoing to uncover the potential of OS to discern the various CLAD phenotypes including bronchiolitis obliterans syndrome (BOS), restrictive allograft syndrome (RAS), and mixed phenotypes⁹.

"OS could be combined with other non-invasive biomarkers, such as computer aided quantitative analysis of chest imaging, or molecular signature of mucosal biopsies to provide a cumulative risk with improved sensitivity and predicative value than individual metrics alone"

-Vasileva et al. 2023

THORASYS

Looking for additional information on oscillometry?
Contact us at Info@thorasys.com

CANADA
6560 de l'Esplanade, Suite 103
Montreal, Quebec H2V 4L5 Canada
1-855-THORASYS | +1-514-384-8555

GERMANY
Tscheulinstraße 21
79331 Teningen, Germany
+49-7641-96-79-353

References: 1. Lundblad, L. K. A. et al. Applications of oscillometry in clinical research and practice. *Can. J. Respir. Crit. Care Sleep Med.* 5, 54–68 (2021). 2. Usmani, O. S. Calling Time on Spirometry: Unlocking the Silent Zone in Acute Rejection after Lung Transplantation. *Am. J. Respir. Crit. Care Med.* 201, 1468–1470 (2020). 3. Cho, E. et al. Airway Oscillometry Detects Spirometric-Silent Episodes of Acute Cellular Rejection. *Am. J. Respir. Crit. Care Med.* 201, 1536–1544 (2020). 4. Deepak, D. et al. Recognition of small airways obstruction in asthma and COPD - The road less travelled. *Journal of Clinical and Diagnostic Research* 11, 5 TE01-TE05 (2017). 5. Van Muylem et al. Role of pulmonary function in the detection of allograft dysfunction after heart-lung transplantation. *Thorax* 52, 643–647 (1997). 6. Benzimra, M. et al. Acute rejection. *Journal of Thoracic Disease* 9, 5440–5457 (2017). 7. Brun, A. L. et al. Lung transplantation: Ct assessment of chronic lung allograft dysfunction (clad). *Diagnostics* vol. 11 (2021). 8. Vasileva, A. et al. Intra-Subject Variance of Respiratory Oscillometry Reflects Graft Injury and is Associated with Acute Rejection and Chronic Lung Allograft Dysfunction (CLAD) Post Lung Transplant (LTx). *J. Hear. Lung Transplant.* 40, S57 (2021). 9. Fu, A. et al. Characterization of chronic lung allograft dysfunction phenotypes using spectral and intrabreath oscillometry. *Front. Physiol.* 13, (2022).